

CD100 antibody

Product Information

Catalog No.:	FNab01407
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	sema domain, immunoglobulin domain(Ig), transmembrane domain(TM) and short cytoplasmic domain,(semaphorin) 4D
Synonyms:	Semaphorin-4D A8 BB18 GR3 SEMA4D C9orf164 CD100 SEMAJ
Observed MW:	150 kDa
Uniprot ID :	Q92854

Application

Reactivity:	Human
Tested Application:	ELISA
Recommended dilution:	None