
anti- CD11a/Integrin alpha-L antibody

Product Information

Catalog No.:	FNab01409
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2D9
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	integrin, alpha L
Synonyms:	CD11A, CD11a/Integrin alpha-L, Integrin alpha L, ITGAL, LFA 1, LFA 1A, LFA1A
Observed MW:	129 kDa
Uniprot ID :	P20701

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:20-1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab01409 (CD11a/Integrin alpha-L antibody) at dilution of 1:400