

ITGAM antibody

Product Information

Catalog No.:	FNab01410
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	integrin, alpha M(complement component 3 receptor 3 subunit)
Synonyms:	Integrin alpha-M CD11 antigen-like family member B CR-3 alpha chain Cell surface glycoprotein MAC-1 subunit alpha Leukocyte adhesion receptor MO1 Neutrophil adherence receptor ITGAM CD11B CR3A
Observed MW:	180 kDa
Uniprot ID :	P11215

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human spleen using FNab01410(CD11B antibody) at dilution of 1:50



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab01410(CD11B/Integrin alpha M Antibody) at dilution of 1:600