

ITGAD antibody

Product Information

Catalog No.:	FNab01414
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	integrin, alpha D
Synonyms:	Integrin alpha-D ADB2 CD11 antigen-like family member D Leukointegrin alpha D ITGAD
Observed MW:	150 kDa
Uniprot ID :	Q13349

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	

Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01414(CD11d antibody) at dilution of 1:800

