

## IL7R antibody

### Product Information

Catalog No.:	FNab01415
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	interleukin 7 receptor
Synonyms:	Interleukin-7 receptor subunit alpha (IL-7 receptor subunit alpha, IL-7R subunit alpha, IL-7R-alpha, IL-7RA) CDw127 IL7R
Observed MW:	56-58 kDa
Uniprot ID :	P16871

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, FC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:2000
Image:	



IP Result of anti-CD127/IL7R (IP:FNab01415, 4ug; Detection:FNab01415 1:1000) with K-562 cells lysate 1200ug.



mouse thymus tissue were subjected to SDS PAGE followed by western blot with FNab01415(CD127/IL7R antibody) at dilution of 1:1500