

SDC1 antibody

Product Information


Catalog No.:	FNab01425
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2G5
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	syndecan 1
Synonyms:	Syndecan-1 (SYND1) SDC1 SDC
Observed MW:	120 kDa
Uniprot ID :	P18827

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsil using FNab01425(SDC1,CD138 antibody) at dilution of 1:50



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab01425(SDC1,CD138 antibody) at dilution of 1:500