

CD40LG antibody

Product Information

Catalog No.:	FNab01432
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	CD40 ligand
Synonyms:	CD40 ligand (CD40-L) T-cell antigen Gp39 TNF-related activation protein (TRAP) Tumor necrosis factor ligand superfamily member 5 CD40 ligand, membrane form CD40 ligand, soluble form (sCD40L) CD40LG CD40L TNFSF5 TRAP
Observed MW:	32 kDa
Uniprot ID :	P29965

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human spleen using FNab01432(CD40LG antibody) at dilution of 1:100



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with FNab01432(CD40LG antibody) at dilution of 1:1000