

CD1A antibody

Product Information

Catalog No.:	FNab01434
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Antigen-presenting protein that binds self and non-self lipid and glycolipid antigens and presents them to T-cell receptors on natural killer T-cells.

Immunogen information

Immunogen:	CD1a molecule
Synonyms:	T-cell surface glycoprotein CD1a T-cell surface antigen T6/Leu-6 (hTa1 thymocyte antigen) CD1A
Observed MW:	65-70 kDa
Uniprot ID :	P06126

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab01434(CD1a antibody) at dilution of 1:800