

CD23 antibody

Product Information


Catalog No.:	FNab01447
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9H0
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	Fc fragment of IgE, low affinity II, receptor for(CD23)
Synonyms:	Low affinity immunoglobulin epsilon Fc receptor BLAST-2 C-type lectin domain family 4 member J Fc-epsilon-RII Immunoglobulin E-binding factor Lymphocyte IgE receptor Low affinity immunoglobulin epsilon Fc receptor membrane-bound form Low affinity immunoglobulin epsilon Fc receptor soluble form FCER2 CD23A CLEC4J FCE2 IGEBF
Observed MW:	38-42 kDa
Uniprot ID :	P06734

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:2000-1:4000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab01447(CD23,FCER2 Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Raji cells were subjected to SDS PAGE followed by western blot with FNab01447(CD23,FCER2 antibody) at dilution of 1:2000