

CD3D antibody

Product Information


Catalog No.:	FNab01457
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8E6
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	CD3d molecule, delta(CD3-TCR complex)
Synonyms:	T-cell surface glycoprotein CD3 delta chain T-cell receptor T3 delta chain CD3D T3D
Observed MW:	24 kDa
Uniprot ID :	P04234

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01457(CD3D Antibody) at dilution of 1:1000