

CD3G antibody

Product Information

Catalog No.:	FNab01460
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	CD3g molecule, gamma(CD3-TCR complex)
Synonyms:	T-cell surface glycoprotein CD3 gamma chain T-cell receptor T3 gamma chain CD3G T3G
Observed MW:	20 kDa
Uniprot ID :	P09693

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB,IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis using FNab01460(CD3G antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01460(CD3G antibody) at dilution of 1:1000