

PECAM1 antibody

Product Information


Catalog No.:	FNab01464
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2D10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	platelet/endothelial cell adhesion molecule
Synonyms:	Platelet endothelial cell adhesion molecule (PECAM-1) EndoCAM GPIIA' PECA1 PECAM1
Observed MW:	130 kDa
Uniprot ID :	P16284

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:2000-1:6000; IHC: 1:500-1:3000
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab01464(CD31 Antibody) at dilution of 1:500



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01464(CD31 antibody) at dilution of 1:1000