

SPN antibody

Product Information


Catalog No.:	FNab01481
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6G11
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	sialophorin
Synonyms:	Leukosialin GPL115 Galactoglycoprotein (GALGP) Leukocyte sialoglycoprotein Sialophorin CD43 cytoplasmic tail (CD43-ct, CD43ct) SPN CD43
Observed MW:	115 kDa
Uniprot ID :	P16150

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab01481 (SPN Antibody) at dilution of 1:100 . Heat mediated antigen retrieved with Tris-EDTA buffer, pH9.0



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01481 (SPN Antibody) at dilution of 1:1000