

CD58 antibody

Product Information

Catalog No.:	FNab01488
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Ligand of the T-lymphocyte CD2 glycoprotein. This interaction is important in mediating thymocyte interactions with thymic epithelial cells, antigen-independent and-dependent interactions of T-lymphocytes with target cells and antigen-presenting cells and the T-lymphocyte rosetting with erythrocytes. In addition, the LFA-3/CD2 interaction may prime response by both the CD2+ and LFA-3+ cells.

Immunogen information

Immunogen:	CD58 molecule
Synonyms:	Lymphocyte function-associated antigen 3 (Ag3) Surface glycoprotein LFA-3 CD58 LFA3
Observed MW:	65-70 kDa
Uniprot ID :	P19256

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Raji cells were subjected to SDS PAGE followed by western blot with FNab01488(CD58 antibody) at dilution of 1:1000