

CD7 antibody

Product Information

Catalog No.:	FNab01495
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Not yet known.

Immunogen information

Immunogen:	CD7 molecule
Synonyms:	T-cell antigen CD7 GP40 T-cell leukemia antigen T-cell surface antigen Leu-9 TP41 CD7
Observed MW:	None
Uniprot ID :	P09564

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab01495(CD7 Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).