

## CD82 antibody

### Product Information

Catalog No.:	FNab01502
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Associates with CD4 or CD8 and delivers costimulatory signals for the TCR/CD3 pathway.

### Immunogen information


Immunogen:	CD82 molecule
Synonyms:	CD82 antigen C33 antigen IA4 Inducible membrane protein R2 Metastasis suppressor Kangai-1 Suppressor of tumorigenicity 6 protein Tetraspanin-27 (Tspan-27) CD82 KAI1 SAR2 ST6 TSPAN27
Observed MW:	30-34 kDa
Uniprot ID :	P27701

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human spleen using FNab01502(CD82 antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with FNab01502(CD82 antibody) at dilution of 1:500