

## CD8B antibody

### Product Information

Catalog No.:	FNab01506
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Identifies cytotoxic/suppressor T-cells that interact with MHC class I bearing targets. CD8 is thought to play a role in the process of T-cell mediated killing.

### Immunogen information

Immunogen:	CD8b molecule
Synonyms:	T-cell surface glycoprotein CD8 beta chain CD8B CD8B1
Observed MW:	37-43 kDa
Uniprot ID :	P10966

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab01506( CD8B Antibody) at dilution of 1:1000