

SLC3A2 antibody

Product Information

Catalog No.:	FNab01511
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	solute carrier family 3(activators of dibasic and neutral amino acid transport), member 2
Synonyms:	Amino acid transporter heavy chain SLC3A2 4F2 cell-surface antigen heavy chain (4F2hc) 4F2 heavy chain antigen Lymphocyte activation antigen 4F2 large subunit Solute carrier family 3 member 2 SLC3A2 MDU1
Observed MW:	85 kDa, 120 kDa
Uniprot ID :	P08195

Application

Reactivity:	Human, rat
Tested Application:	ELISA, WB, IHC, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab01511(SLC3A2 Antibody) at dilution of 1:500



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01511(SLC3A2 antibody) at dilution of 1:500