

## CDK19 antibody

### Product Information

Catalog No.:	FNab01526
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	cell division cycle 2-like 6(CDK8-like)
Synonyms:	Cyclin-dependent kinase 19 CDC2-related protein kinase 6 Cell division cycle 2-like protein kinase 6 Cell division protein kinase 19 Cyclin-dependent kinase 11 Death-preventing kinase CDK19 CDC2L6 CDK11 KIAA1028
Observed MW:	60-64 kDa
Uniprot ID :	Q9BWU1

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse eye tissue were subjected to SDS PAGE followed by western blot with FNab01526(CDC2L6 antibody) at dilution of 1:300