

## CDC42EP2 antibody

### Product Information

Catalog No.:	FNab01531
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Probably involved in the organization of the actin cytoskeleton. May act downstream of CDC42 to induce actin filament assembly leading to cell shape changes. Induces pseudopodia formation in fibroblasts in a CDC42-dependent manner.

### Immunogen information

Immunogen:	CDC42 effector protein(Rho GTPase binding) 2
Synonyms:	Cdc42 effector protein 2 Binder of Rho GTPases 1 CDC42EP2 BORG1 CEP2
Observed MW:	22 kDa
Uniprot ID :	O14613

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab01531(CDC42EP2 antibody) at dilution of 1:200