

CDC42EP4 antibody

Product Information

Catalog No.:	FNab01532
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probably involved in the organization of the actin cytoskeleton. May act downstream of CDC42 to induce actin filament assembly leading to cell shape changes. Induces pseudopodia formation, when overexpressed in fibroblasts.

Immunogen information

Immunogen:	CDC42 effector protein(Rho GTPase binding) 4
Synonyms:	Cdc42 effector protein 4 Binder of Rho GTPases 4 CDC42EP4 BORG4 CEP4
Observed MW:	43-45 kDa
Uniprot ID :	Q9H3Q1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with 1FNab01532(CDC42EP4 antibody) at dilution of 1:1000