

CDC6 antibody

Product Information

Catalog No.:	FNab01536
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the initiation of DNA replication. Also participates in checkpoint controls that ensure DNA replication is completed before mitosis is initiated.

Immunogen information


Immunogen:	cell division cycle 6 homolog
Synonyms:	Cell division control protein 6 homolog CDC6-related protein Cdc18-related protein (HsCdc18) p62(cdc6) (HsCDC6) CDC6 CDC18L
Observed MW:	63 kDa
Uniprot ID :	Q99741

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab01536(CDC6 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01536(CDC6 antibody) at dilution of 1:400