

CDC7 antibody

Product Information

Catalog No.:	FNab01538
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	cell division cycle 7 homolog
Synonyms:	Cell division cycle 7-related protein kinase (CDC7-related kinase, HsCdc7, huCdc7) CDC7 CDC7L1
Observed MW:	64 kDa
Uniprot ID :	O00311

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:10-1:100
Image:	



H-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab01538(Cdc7 Kinase antibody) at dilution of 1:500