

CDK2AP1 antibody

Product Information

Catalog No.:	FNab01555
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


specific inhibitor of the cell-cycle kinase CDK2.

Immunogen information

Immunogen:	cyclin-dependent kinase 2 associated protein 1
Synonyms:	Cyclin-dependent kinase 2-associated protein 1 (CDK2-associated protein 1) Deleted in oral cancer 1 (DOC-1) Putative oral cancer suppressor CDK2AP1 CDKAP1 DOC1
Observed MW:	12 kDa
Uniprot ID :	O14519

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab01555(CDK2AP1 antibody) at dilution of 1:500