

CDK5RAP1 antibody

Product Information

Catalog No.:	FNab01559
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Specifically inhibits CDK5 activation by CDK5R1.

Immunogen information


Immunogen:	CDK5 regulatory subunit associated protein 1
Synonyms:	Mitochondrial tRNA methylthiotransferase CDK5RAP1 CDK5 activator-binding protein C42 CDK5 regulatory subunit-associated protein 1 mt-tRNA-2-methylthio-N6-dimethylallyl-adenosine synthase mt-tRNA-N6-(dimethylallyl)adenosine(37) methylthiotransferase CDK5RAP1 C20orf34
Observed MW:	68 kDa
Uniprot ID :	Q96SZ6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas using FNab01559(CDK5RAP1 antibody) at dilution of 1:100



human heart tissue were subjected to SDS PAGE followed by western blot with FNab01559(CDK5RAP1 antibody) at dilution of 1:2000