

## CEACAM21 antibody

### Product Information

Catalog No.:	FNab01577
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	carcinoembryonic antigen-related cell adhesion molecule 21
Synonyms:	Carcinoembryonic antigen-related cell adhesion molecule 21 CEACAM21
Observed MW:	32, 22 kDa
Uniprot ID :	Q3KPI0

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01577(CEACAM21 antibody) at dilution of 1:500