

CENPJ antibody

Product Information

Catalog No.:	FNab01587
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays an important role in cell division and centrosome function by participating in centriole duplication. Inhibits microtubule nucleation from the centrosome. Required for centriole elongation and for STIL-mediated centriole amplification(PubMed:22020124).

Immunogen information


Immunogen:	centromere protein J
Synonyms:	Centromere protein J (CENP-J) Centrosomal P4.1-associated protein LAG-3-associated protein LYST-interacting protein 1 CENPJ CPAP LAP LIP1
Observed MW:	153 kDa
Uniprot ID :	Q9HC77

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IF: 1:10-1:100
Image:	



Immunofluorescent analysis of HepG2 cells using FNab01587(CENPJ Antibody) at dilution of 1:30 and Rhodamine-Goat anti-Rabbit IgG



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab01587(CENPJ antibody) at dilution of 1:500