

CEP170P1 antibody

Product Information

Catalog No.:	FNab01599
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	centrosomal protein 170kDa-like
Synonyms:	Cep170-like protein CEP170 pseudogene 1 CEP170P1 CEP170L KIAA0470L
Observed MW:	None
Uniprot ID :	Q96L14

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	

HeLa cells were subjected to SDS PAGE followed by western blot with FNab01599(CEP170L) antibody at dilution of 1:1000

