

CEP55 antibody

Product Information

Catalog No.:	FNab01605
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role in mitotic exit and cytokinesis. Not required for microtubule nucleation. Recruits PD CD61P and TSG101 to midbody during cytokinesis.

Immunogen information


Immunogen:	centrosomal protein 55kDa
Synonyms:	Centrosomal protein of 55 kDa (Cep55) Up-regulated in colon cancer 6 CEP55 C10orf3 URCC6
Observed MW:	47-55 kDa
Uniprot ID :	Q53EZ4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab01605(CEP55 Antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab01605(CEP55 Antibody) at dilution of 1:1000