

CEP70 antibody

Product Information

Catalog No.:	FNab01608
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role in the organization of both preexisting and nascent microtubules in interphase cells. During mitosis, required for the organization and orientation of the mitotic spindle.

Immunogen information


Immunogen:	centrosomal protein 70kDa
Synonyms:	Centrosomal protein of 70 kDa (Cep70) p10-binding protein CEP70 BITE
Observed MW:	25-28 kDa
Uniprot ID :	Q8NHQ1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary using FNab01608(CEP70 antibody) at dilution of 1:50



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab01608(CEP70 antibody) at dilution of 1:600