

CERCAM antibody

Product Information

Catalog No.:	FNab01612
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probable cell adhesion protein involved in leukocyte transmigration across the blood-brain barrier. Has apparently no beta-galactosyltransferase activity.

Immunogen information

Immunogen:	cerebral endothelial cell adhesion molecule
Synonyms:	Inactive glycosyltransferase 25 family member 3 Cerebral endothelial cell adhesion molecule CERCAM CEECAM1 GLT25D3 KIAA1502
Observed MW:	68 kDa, 48 kDa
Uniprot ID :	Q5T4B2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab01612(CERCAM antibody) at dilution of 1:1500