

CP antibody

Product Information


Catalog No.:	FNab01614
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8F3
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	ceruloplasmin(ferroxidase)
Synonyms:	Ceruloplasmin Ferroxidase CP
Observed MW:	140 kDa
Uniprot ID :	P00450

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using FNab01614(CP Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab01614(CP Antibody) at dilution of 1:1000