

CES5A-Specific antibody

Product Information

Catalog No.:	FNab01618
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	carboxylesterase 7
Synonyms:	Carboxylesterase 5A Carboxylesterase-like urinary excreted protein homolog (Cauxin) CES5A CES7
Observed MW:	58 kDa
Uniprot ID :	Q6NT32

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:1000; IF: 1:20-1:200
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab01618(CES7-Specific antibody) at dilution of 1:300