

CFLAR antibody

Product Information

Catalog No.:	FNab01621
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	CASP8 and FADD-like apoptosis regulator
Synonyms:	CASP8 and FADD-like apoptosis regulator Caspase homolog (CASH) Caspase-eight-related protein (Casper) Caspase-like apoptosis regulatory protein (CLARP) Cellular FLICE-like inhibitory protein (c-FLIP) FADD-like antiapoptotic molecule 1 (FLAME-1) CASH CASP8AP1 CLARP MRIT
Observed MW:	30-35 kDa, 50-55 kDa
Uniprot ID :	O15519

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle slide using FNab01621(CFLAR Antibody) at dilution of 1:50



IP Result of anti-CFLAR/FLIP (IP:FNab01621, 4ug; Detection:FNab01621 1:500) with Jurkat cells lysate 4000ug.



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01621(CFLAR antibody) at dilution of 1:1000