

CHD2 antibody

Product Information

Catalog No.:	FNab01641
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


DNA-binding helicase that specifically binds to the promoter of target genes, leading to chromatin remodeling, possibly by promoting deposition of histone H3.3. Involved in myogenesis via interaction with MYOD1: binds to myogenic gene regulatory sequences and mediates incorporation of histone H3.3 prior to the onset of myogenic gene expression, promoting their expression(By similarity).

Immunogen information


Immunogen:	chromodomain helicase DNA binding protein 2
Synonyms:	Chromodomain-helicase-DNA-binding protein 2 (CHD-2) ATP-dependent helicase CHD2 CHD2
Observed MW:	95 kDa
Uniprot ID :	O14647

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab01641(CHD2 antibody) at dilution of 1:50



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab01641(CHD2 antibody) at dilution of 1:500