

CHD4 antibody

Product Information

Catalog No.:	FNab01643
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones.

Immunogen information


Immunogen:	chromodomain helicase DNA binding protein 4
Synonyms:	Chromodomain-helicase-DNA-binding protein 4 (CHD-4) ATP-dependent helicase CHD4 Mi-2 autoantigen 218 kDa protein Mi2-beta CHD4
Observed MW:	260 kDa
Uniprot ID :	Q14839

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cervical cancer using FNab01643(CHD4 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01643(CHD4 antibody) at dilution of 1:500