

## CHD4 antibody

### Product Information

Catalog No.:	FNab01644
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	7E0
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones.

### Immunogen information

Immunogen:	chromodomain helicase DNA binding protein 4
Synonyms:	Chromodomain-helicase-DNA-binding protein 4 (CHD-4) ATP-dependent helicase CHD4 Mi-2 autoantigen 218 kDa protein Mi2-beta CHD4
Observed MW:	260 kDa
Uniprot ID :	Q14839

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab01644( CHD4 Antibody) at dilution of 1:1000