

## CHRNA5 antibody

### Product Information

Catalog No.:	FNab01677
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


After binding acetylcholine, the AChR responds by an extensive change in conformation that affects all subunits and leads to opening of an ion-conducting channel across the plasma membrane.

### Immunogen information

Immunogen:	cholinergic receptor, nicotinic, alpha 5
Synonyms:	Neuronal acetylcholine receptor subunit alpha-5 CHRNA5 NACHRA5
Observed MW:	53 kDa
Uniprot ID :	P30532

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01677(CHRNA5 antibody) at dilution of 1:500