

CHRNA4 antibody

Product Information

Catalog No.:	FNab01684
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


After binding acetylcholine, the AChR responds by an extensive change in conformation that affects all subunits and leads to opening of an ion-conducting channel across the plasma membrane.

Immunogen information

Immunogen:	cholinergic receptor, nicotinic, beta 4
Synonyms:	Neuronal acetylcholine receptor subunit beta-4 CHRNA4
Observed MW:	48-60 kDa
Uniprot ID :	P30926

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab01684(CHRNB4 antibody) at dilution of 1:500