

## SLC5A7 antibody

### Product Information

Catalog No.:	FNab01700
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	solute carrier family 5(choline transporter), member 7
Synonyms:	High affinity choline transporter 1 (hCHT1) Hemicholinium-3-sensitive choline transporter (CHT) Solute carrier family 5 member 7 SLC5A7 CHT1
Observed MW:	65-70 kDa
Uniprot ID :	Q9GZV3

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain slide using FNab01700(SLC5A7 Antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01700(SLC5A7 antibody) at dilution of 1:300