

CILP2 antibody

Product Information

Catalog No.:	FNab01711
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in cartilage scaffolding.

Immunogen information

Immunogen:	cartilage intermediate layer protein 2
Synonyms:	Cartilage intermediate layer protein 2 (CILP-2) Cartilage intermediate layer protein 2 C1 Cartilage intermediate layer protein 2 C2 CILP2
Observed MW:	90 kDa
Uniprot ID :	Q8IUL8

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab01711(CILP2 antibody) at dilution of 1:500