

SARNP antibody

Product Information

Catalog No.:	FNab01714
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	cytokine induced protein 29 kDa
Synonyms:	SAP domain-containing ribonucleoprotein Cytokine-induced protein of 29 kDa Nuclear protein Hcc-1 Proliferation-associated cytokine-inducible protein CIP29 SARNP HCC1
Observed MW:	29-32 kDa
Uniprot ID :	P82979

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab01714(CIP29 Antibody) at dilution of 1:600