

Utp4 antibody

Product Information

Catalog No.:	FNab01716
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	cirrhosis, autosomal recessive 1A(cirhin)
Synonyms:	Utp4 Cirh1a
Observed MW:	70-77 kDa
Uniprot ID :	Q5I0C9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab01716(Cirhin antibody) at dilution of 1:500