

CKS2 antibody

Product Information

Catalog No.:	FNab01729
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function.

Immunogen information

Immunogen:	CDC28 protein kinase regulatory subunit 2
Synonyms:	Cyclin-dependent kinases regulatory subunit 2 (CKS-2) CKS2
Observed MW:	10 kDa
Uniprot ID :	P33552

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Rat live tissue were subjected to SDS PAGE followed by western blot with FNab01729(CKS2 antibody) at dilution of 1:1000