

## CLCN1 antibody

### Product Information

Catalog No.:	FNab01743
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Voltage-gated chloride channel. Chloride channels have several functions including the regulation of cell volume; membrane potential stabilization, signal transduction and transepithelial transport.

### Immunogen information

Immunogen:	chloride channel 1, skeletal muscle
Synonyms:	Chloride channel protein 1 (ClC-1) Chloride channel protein, skeletal muscle CLCN1 CLC1
Observed MW:	109 kDa
Uniprot ID :	P35523

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



293T cells were subjected to SDS PAGE followed by western blot with FNab01743(CLCN1 antibody) at dilution of 1:1000