

CLOCK antibody

Product Information


Catalog No.:	FNab01769
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	clock homolog(Mouse)
Synonyms:	Circadian locomoter output cycles protein kaput (hCLOCK) Class E basic helix-loop-helix protein 8 (bHLHe8) CLOCK BHLHE8 KIAA0334
Observed MW:	110 kDa
Uniprot ID :	O15516

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP, IF
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



IP result of anti-CLOCK(IP:FNab01769, 4ug;
Detection:FNab01769 1:300) with HeLa cells
lysate 1600 ug



Various lysates were subjected to SDS PAGE
followed by western blot with
FNab01769(CLOCK antibody) at dilution of 1:500