

CLPS antibody

Product Information

Catalog No.:	FNab01773
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Colipase is a cofactor of pancreatic lipase. It allows the lipase to anchor itself to the lipid-water interface. Without colipase the enzyme is washed off by bile salts, which have an inhibitory effect on the lipase. Enterostatin has a biological activity as a satiety signal.

Immunogen information

Immunogen:	colipase, pancreatic
Synonyms:	Colipase CLPS
Observed MW:	12 kDa
Uniprot ID :	P04118

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with FNab01773(CLPS antibody) at dilution of 1:600