

CLUAP1 antibody

Product Information

Catalog No.:	FNab01780
Size:	100μg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for cilia biogenesis. Appears to function within the multiple intraflagellar transport complex B(IFT-B). Key regulator of hedgehog signaling(By similarity).

Immunogen information

Immunogen:	clusterin associated protein 1
Synonyms:	Clusterin-associated protein 1 Qilin CLUAP1 KIAA0643
Observed MW:	53 kDa
Uniprot ID :	Q96AJ1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab01780(CLUAP1 antibody) at dilution of 1:600