

## CNKSR1 antibody

### Product Information

Catalog No.:	FNab01798
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May function as an adapter protein or regulator of Ras signaling pathways.

### Immunogen information


Immunogen:	connector enhancer of kinase suppressor of Ras 1
Synonyms:	Connector enhancer of kinase suppressor of ras 1 (Connector enhancer of KSR 1) CNK homolog protein 1 (CNK1, hCNK1) Connector enhancer of KSR-like CNKSR1 CNK1
Observed MW:	90-100 kDa
Uniprot ID :	Q969H4

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using FNab01798(CNKSR1 Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



HT-29 cells were subjected to SDS PAGE followed by western blot with FNab01798(CNKSR1 Antibody) at dilution of 1:600