

## CNOT4 antibody

### Product Information

Catalog No.:	FNab01809
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Has E3 ubiquitin ligase activity. Involved in activation of the JAK/STAT pathway.

### Immunogen information


Immunogen:	CCR4-NOT transcription complex, subunit 4
Synonyms:	CCR4-NOT transcription complex subunit 4 CCR4-associated factor 4 E3 ubiquitin-protein ligase CNOT4 Potential transcriptional repressor NOT4Hp RING-type E3 ubiquitin transferase CNOT4 CNOT4 NOT4
Observed MW:	64 kDa
Uniprot ID :	O95628

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab01809(CNOT4 antibody) at dilution of 1:50



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab01809(CNOT4 antibody) at dilution of 1:500